

OPERATORY SEATING

 Belmont®

Model 100
Doctor's
Memory Foam Stool



Model 101
Assistant's
Memory Foam Stool



Model 100

Doctor's

Memory Foam Stool

- Memory Foam Core Seat Cushion
- Individual Levers For Seat Height & Backrest Tilt
- Adjustable Seat Height Range: 18"-23"
- Vertical Backrest Height Adjustment
- Smooth, Plated 5-Spoke Base
- Metal-Hooded Casters With Dual Black Wheels

Upholstery:

Seamed Style

Duration (Belmont Standard Colors)

Options:

Ultraleather™ Upholstery

Short Cylinder for 16"-21" Seat Height Range

Model 101

Assistant's

Memory Foam Stool

- Memory Foam Core Seat Cushion
- Horizontal Backrest Adjustment With 0 Degree Tilt
- Height-Adjustable, Ratcheting Armrest Support
- Adjustable Seat Height Range: 22"-28"
- Smooth, Plated 5-Spoke Base
- Metal-Hooded Casters With Dual Black Wheels
- Height-Adjustable 18" Diameter Tubular Foot Ring

Note:

Specify Armrest Support for Right-Handed Or Left-Handed User When Ordering

Upholstery:



Seamed Style

Duration (Belmont Standard Colors)

Options:

Ultraleather™ Upholstery

Long Cylinder for 24"-33" Seat Height Range



Luxurious Memory Foam Seat Cushions Conform To Your Body Shape
Seat Cushion Shape Recovers When No Weight Is Applied

Belmont®

Belmont Equipment
A Division of Takara Belmont USA, Inc.
101 Belmont Drive
Somerset, NJ 08873
(732) 469-5000

(800)223-1192
www.belmontequip.com