

# OPERATORY SEATING

 Belmont®



# DENTIST SEATING

**090**

Doctor's  
Stool



- Ergonomic chair control with a single lock lever for seat and backrest tilt adjustments
- Satin brushed aluminum finish on top of base with painted black spokes
- Chrome hooded casters with dual black wheels
- Seamless upholstery
- Backrest horizontal adjustment
- Degree tilt
- 18" Diameter tubular footring
- Ratcheting body support mechanism
- Right and left hand models (Model 091)

Model	Seat Height Range	Seat Tilt
090	18"– 24" (Standard)	Yes
	16.5"– 21" (Optional)	

Model	Seat Height Range	Seat Tilt
091	20"– 26" (Standard)	No
	22"– 31" (Optional)	

**091**

Assistant's  
Stool



## Upholstery

Standard	Morbern Duration Vinyl
Ultraleather™	Ultrafabrics, LLC polyurethane at an upcharge

**Belmont®**

Belmont Equipment  
A Division of Takara Belmont USA, Inc.  
101 Belmont Drive,  
Somerset, New Jersey 08873  
(732) 469-5000

**1-800-223-1192**  
**www.belmontequip.com**

Copyright © 2001 Belmont Equipment  
BRH00008 REV103103