

E-120 Quolis[®] STERILIZATION CENTERS

E-Series Cabinetry

Belmont
www.belmontequip.com

E-120ECS Unit Shown in Folkstone / Sheffield Cambria Top
Overall Dimensions: 121" (L) x 84" (H) x 27-3/8" (D)



High Quality Custom Designs to Fit Your Practice's Needs!

E-120 Quolis® STERILIZATION CENTERS

E-Series Cabinetry

St. Louis, MO: 1-888-833-3682
Somerset, NJ: 1-800-223-1192

Belmont
www.belmontequip.com

Unsterile Upper Instrument Cabinet w/ Glass Shelves and Red Light

Upper Storage w/ (2) Glove Dispenser

Upper Storage

Autoclave Sterilizer Tower w/ Upper Storage, (2) Pull-Out Shelves, (2) 3" Wrap Drawers, (2) Duplex GFI Outlets, Lower Storage, & Slide-Out Platform

FEATURES:

- Soft-Closing Mechanism on Doors and Drawers*
- Stainless Steel Pulls on All Doors and Drawers
- 3 mm PVC Edging on All Doors and Drawers
- Snap-On Concealed European Hinges on All Doors

Standard E-120 Overall Dimensions: 119-1/4" (L) x 84" (H) x 27-3/8" (D)
* 3" Wrap Drawers in Autoclave Tower Do Not Have this Feature



Pull-Out Rotating Soft-Closing Instrument Tower w/Light



Pull-Out Waste Container




Instrument Cabinet w/ Glass Shelves and Light



Autoclave Tower Slide-Out Platform

OPTIONS:

- Cambria™ Quartz Countertop
- Task Lighting
- White Faucet and Soap Dispenser
- Undermount Stainless or White Sink
- Sharps Disposal Container
- Ultrasonic Cut-Out in Countertop
- Finished Back and End Panel(s)



Optional Telescopic Faucet

E-120 Unit Shown in Folkstone / Sheffield Cambria Top

Cup and Towel Dispenser w/ Optional Task Light

Chrome Faucet w/ Stainless Steel Sink & Soap Dispenser

Trash Chute & Optional Sharps Container

Lower Storage Cabinet for Ultrasonic Unit

Pull-Out Waste Container and Full Extension Soft-Close Tandembox Drawer

Lower Sink Storage with Plumbing Access

Foot Controlled Tapmaster™ Sink Activation

(2) 5", (1) 7", & (1) 10" Full Extension Soft-Close Tandembox Drawers


Sanitized Pull-Out Rotating Soft-Closing Instrument Tower w/ Blue Light



Visualize



Design



Develop



Let Belmont Equipment Assist You In Designing the Perfect Sterilization Center for your Practice!

Items may be shown with optional equipment. Please consult Customer Service for complete details. Belmont © reserves the right to change standards, options, and prices without notice.

©2009 Belmont Equipment

BRH00045 ver 1.0