

# PRO II

## Dental Chair

**Consultation, Case Presentation,  
Prosthetics and Handicapped Accessible**

**The latest generation of PRO II chairs with**

- ❑ Articulating dual axis headrest
- ❑ Surfaces are easily accessible to wipe-down with disinfectant
- ❑ Protection from disease transmission by finger to switch contamination
- ❑ Eye to oral cavity contact can be maintained during the procedure
- ❑ Communication between you and your patient is accomplished in a relaxed and professional manner; your patient is at ease
- ❑ 12 positioning programs to meet your individual requirements



**Model 037-S**

# PRO II Dental Chair

Model **037-S**

Protection For Professionals

### Articulating Dual Axis Prosthetic headrest

Adjusts to match patient's anatomy. Locks securely in any desired position.

### The Wing-Design Upholstery

Insures patient shoulder support.

### Armrest Swing-Out

For ease in patient entry and firm support on exiting the chair.

### Choice of Upholstery

Available with Duration™ or Ultraleather™ at no additional charge.



### Thin Backrest

for maximum body positioning and access to the oral cavity.

### The Lumbar Support

is anatomically correct in all patient positions.

### Roll-Out Legrest

The legrest extends when the chair is placed in a pre-set position.

### Program Panel

The dentist can set all twelve pre-set positions by "dialing-in" the location for the seat and back.

### Foot Controls

Provide "hands-free" operation for selection of twelve pre-set positions, patient dismissal, and back/base movements.



### The Rotation is 110° L&R

For face-to-face patient consultation.

### 5 Year Warranty

On mechanical and electrical parts; it applies when a Belmont chair is purchased directly from an authorized dealer; and, is installed and serviced by an employee of an authorized dealer.

The purchase of the PRO II chair may qualify you for the disabled access tax credit (IRC44). Your tax advisor should be consulted to determine if the credit is applicable to your purchase.

Belmont reserves the right to change specifications and pricing without notice.  
All rights reserved: ©Copyright Belmont Equipment.  
BRH00013 REV120103

 **Belmont**®

Belmont Equipment  
A Division of Takara Belmont USA, Inc.  
101 Belmont Drive, Somerset, NJ 08873  
(732) 469-5000

**1-800-223-1192**  
[www.belmontequip.com](http://www.belmontequip.com)