

V-0607-060,060

090-16, 090-U-16, 090-V-16




WP/N = SUPPLIED AS COMPONENT OF
PART NUMBER LISTED WITHOUT ASTERISK

21	PLATE KNOB	97781	1	97881	1
20	HARDWARE KIT	97657	1	97657	1
19	UPRIGHT, BR	50939	1	50939	1
17	FLAT WASHER	07618	2	07618	2
16	SPRING	10-310	1	10-310	1
15	SCREW, MACHINE	98014	2	98014	2
14	SCREW, WOOD	98015	2	98015	2
13	BACKREST HINGE	10-311	1	10-311	1
12	SHROUD, BACKREST	98013	1	98013	1
11	SERVICE PART UPH BACK w/SHROUD	98018	1	98018	1
10	MACHINE SCREW	07489	8	07489	8
9	ACTUATOR PADDLE	97781	1	97781	1
8	ACTUATOR LEVER	97781	1	97781	1
7	E-RING	97781	1	97781	1
6	CONTROL PLATE	97781	1	97781	1
5	SHROUD, SEAT	98012	1	98012	1
4	SERVICE PART UPH SEAT w/SHROUD	98017	1	98017	1
3	CYLINDER, W/COLUMN	08994	1	98126	1
2	BASE	98011	1	98011	1
1	CASTER	10-044	5	10-044	5
ITEM NO	BQM. ITEM DESCRIPTION	PART#	QTY	PART#	QTY
	DOCTOR STOOL	090-CU.VY		090-CU.VY-16	
	MODEL NUMBER				

ALL BELMONT ASST'S STOOLS:

091, 091-22, 091R, AND 019R-22



*P/N = SUPPLIED AS COMPONENT OF
PART NUMBER LISTED WITHOUT ASTERISK

44	KNOB	#98010	1	1	1	1	1
43	WASHER	#98010	1	1	1	1	1
42	BACKREST SHROUD	98013	1	1	1	1	1
41	UPRIGHT, BR	09328	1	1	1	1	1
40	ADJ. INT. BR. MCA #98010	98018	1	1	1	1	1
39	MACHINE SCREW	98015	7	7	7	7	7
38	MACHINE SCREW	98014	2	2	2	2	2
35	SPRING, BR HINGE	10-310	1	1	1	1	1
34	BRACKET, BR HINGE	10-311	1	1	1	1	1
29	FLAT WASHER	07618	4	4	4	4	4
28	MACHINE SCREW	07489	8	8	8	8	8
26	ACTUATOR PADDLE	#98010	1	1	1	1	1
25	ACTUATOR LEVER	#98010	1	1	1	1	1
24	E-RING	#98010	1	1	1	1	1
23	CONTROL PLATE	98010	1	1	1	1	1
22	SEAT SHROUD	98022	1	1	1	1	1
21	SEAT INT. BR. SHROUD #98010	98029	1	1	1	1	1
19	SHROUD	97661	1	1	1	1	1
18	MACHINE SCREW	97633	3	3	3	3	3
17	PATCHET, LH	97689	1	1	1	1	1
17	PATCHET, RH	97693	1	1	1	1	1
16	SHROUD INSERT	97661-1	1	1	1	1	1
12	SEAT INT. BR. SHROUD #98010	98019	1	1	1	1	1
10	FOOTRING	97672	1	1	1	1	1
9	SWINGARM ASSY	97682	1	1	1	1	1
8	KEY NUT	758094	1	1	1	1	1
7	KNOB	08083	1	1	1	1	1
3	CYLINDER, V/COLUMN	50945	1	1	1	1	1
3	CYLINDER, V/COLUMN	98127	1	1	1	1	1
2	BASE	98011	1	1	1	1	1
1	CASTER	10-044	5	5	5	5	5
MODEL NUMBER							
P/N							
QTY							
QTY							
QTY							
QTY							