




80B & 81B Dental Stools



■ 80B Doctor Stool

Height Range: 18"-24"



Clean Form and Function

Coordinated seat and backrest tilt permits the operator to maintain an efficient and comfortable working posture.

Backrest Height

Backrest cushion height adjustment provides optimal lumbar support to suit the preference of each operator.

Smooth Movement

All stools are equipped with hard-floor surface type casters to provide smooth movement with chrome metal hoods for protection against accidental damage.

■ 81B Assistant Stool

Height Range: 20"-26"

Body support armrest is available for left- or right-handed assistants.



Adjustable Backrest

Backrest in-and-out adjustment promotes a good working posture by positioning the assistant upright and forward, towards the body support cushion.

Footring Capabilities

To better establish a comfortable and efficient working posture for assistants, the height adjustable tubular footring is used to relieve pressure from the back of the legs.

Body Support

Includes a height adjustable ratcheting body support to provide comfort and relieve stress to arms and back during procedures.

Belmont Equipment
101 Belmont Drive
Somerset, NJ 08873
(800) 223-1192
www.belmontequip.com

Takara Company Canada Ltd.
2076 South Sheridan Way
Mississauga, ON
L5J 2M4, Canada
(800) 268-5351
dental.takarabelmont.ca

 **Belmont**
Delivering Value in Quality

@BelmontEquipUS

