



NOTES:

- 1.) LOCAL REGULATIONS PROVIDE THAT LICENSED PLUMBERS AND ELECTRICIANS SHALL INSTALL UTILITIES
 - 2.) MAKE SURE ALL PLUMBING CONFORMS TO PREVAILING LOCAL CODES
 - 3.) PLUMBING CENTRAL VACUUM TO THE UTILITY CENTER SHOULD BE SPECIFIED BY THE CENTRAL VACUUM SUPPLIER AND TERMINATED IN UTILITY CENTER WITH 3/4" O.D. TUBE PERPENDICULAR TO THE FLOOR, SIMILAR TO THE DRAIN CONNECTION
 - 4.) GRAVITY DRAIN HAS 3/4" NOMINAL PIPE PROTRUDING 1." FROM FINISHED FLOOR
- NOTE: PLACE TRAP IN LINE TO CONFORM WITH LOCAL CODES. SUPPLIED BY CONTRACTOR

ECO-Sys Delivery Systems

Junction Box Template



101 Belmont Drive, Somerset, NJ 08873
 Phone: (800)223-1192 FAX: (732)35-1035
 www.belmontequip.com