



Top of X-ray Horizontal Arm Extends Approximately 3-1/4" Above the Top Edge of the Wall Plate

17-13/16" (REF.)

11-11/16"

5-13/16"

Ø1" Hole for 4 conductor cable

Ø1" Hole for AC 120V

NOTE: Level top and side of plate prior to securing to the wall

Drill (4) Ø 3/16" pilot holes

Lower Edge of Plastic Cover

Lower Edge of Plastic Cover

Minimum Mounting Height from Floor to Bottom Edge of Plate: 34-1/2"

FLOOR

**Bel-Ray II (Model 097)
PHOT-xlIs (Model 505)**
Wall Mount Template
Belmont
1-800-223-1192