




090 & 091 Dental Stools



■ 090 Doctor Stool

Standard seat height ranges from 18"-24"
(Optional height 16.5"-21")



Clean Form and Function

Lever control for coordinated seat/backrest tilt adjustment optimizes comfort and helps align posture.

Backrest Height

Backrest cushion height adjustment provides optimal lumbar support to suit the preference of each operator during procedures and consultations.

Smooth Movement

All stools are equipped with hard-floor surface type casters to provide smooth movement with chrome metal hoods for protection against accidental damage.

■ 091 Assistant Stool

Standard seat height ranges from 20"-26"
(optional height 22"-31")

Body support armrest is available for left- or right-handed assistants.



Adjustable Backrest

Backrest in-and-out adjustment promotes a good working posture by positioning the assistant upright and forward, towards the body support cushion.

Footring Capabilities

To better establish a comfortable and efficient working posture, the height adjustable tubular footring is used to accommodate Assistants that are seated in an elevated working position by providing a footrest that relieves pressure from the back of the legs.

Body Support

Ratcheting body support armrest with adjustable height cushion provides comfort and relieves stress to the arms and back during procedures.

Belmont Equipment
101 Belmont Drive
Somerset, NJ 08873
(800) 223-1192
www.belmontequip.com

Takara Company Canada Ltd.
2076 South Sheridan Way
Mississauga, ON
L5J 2M4, Canada
(800) 268-5351
dental.takarabelmont.ca

Belmont
Delivering Value in Quality

@BelmontEquipUS

