



BDU

Delivery System



BDU Delivery System

A true example of reliability, flexibility, and durability.



BDU570
Swing Mount



BDU500
Over The Patient

Easy Placement

Balance arm with pneumatic brake places the delivery head where it's needed and locks it in place.

Operatory Footprint

Small delivery head takes up minimal space with a flat top to secure a stainless steel tray and non-slip pad to lay out instrumentation.

Integrated Features

Pre-wired to support unit mounted touchpad and other add-on accessories.

Versatility

Mounts to Belmont's B50N and 037N dental chairs, also works with Bel-Halo, Clesta LED and Clesta Halogen Lights.

Water System

City/bottle water system with water bottle quick disconnect attached to PMU or Swing Arm.

Belmont Equipment
101 Belmont Drive
Somerset, NJ 08873
(800) 223-1192
www.belmontequip.com

Takara Company Canada Ltd.
2076 South Sheridan Way
Mississauga, ON
L5J 2M4, Canada
(800) 268-5351
dental.takarabelmont.ca

 **Belmont**
Delivering Value in Quality

@BelmontEquipUS

