



PRO III

Versatile Dental
Chair Package

PRO III
Dental Chair

Evogue
Delivery System



Index

Pro III

Versatile Dental Chair Package

01 **Pro III** Dental Chair

02 **Evogue** Delivery System

01

Pro III
Dental Chair

A Chair that Supports You from Every Angle.

**Experience the ultimate in dental
comfort and functionality with the
Pro III dental chair.**

Introducing the Pro III dental chair—a game-changer in innovation, functionality, and comfort. Tailored for both practitioners and patients, it merges advanced ergonomics with user-friendly features, setting a new benchmark in dental care. The Pro III offers superior support, ensuring comfort and confidence at every appointment.



Folding Legrest & Seat Rotation

The Pro III chair elevates the dental experience with features that set it apart from the rest. Its easy patient entry, space-saving efficiency, and unique design elements, such as swing-out armrests, automatic legrest extension, and a dual-articulating headrest, surpass standard offerings, providing unmatched comfort for all patients.



Upright Position

Reclined Position



Upright Position

The extendable legrest and folding armrests facilitate easy patient entry before smoothly transitioning to a traditional dental setup for treatment.



Reclined Position

Built-in foot control levers reduce floor clutter around the chair and allow operators to transition seated patients to the desired treatment position.



Exceptional Capabilities, All-Around Versatility

The Pro III chair offers unparalleled flexibility with its 210-degree rotation range, allowing effortless positioning for easy patient entry, relaxed consultations, and seamless examinations from multiple angles.





10:00 Position



Doctors and staff can adjust the chair to create a welcoming environment and facilitate face-to-face consultations, reducing patient anxiety.



9:00 Position



Doctors can pivot patients to an optimal position, enhancing interactions and fostering a comfortable trusting environment.



6:00 Position



The unique folding legrest allows doctors to access patients seated upright directly from the front, promoting a natural, relaxed patient position for effective jaw alignment assessment and restorative procedures.


Reassuring Qualities

The Pro III chair is designed for optimal accessibility, with rotating armrests that ease patient entry and exit, reducing strain for both patients and staff. Built for durability and stability, the chair supports patients of all ages and mobility levels, with upright capabilities and rotating armrests providing easy, unobstructed access.




Oral Cavity Access

The chair's dual-articulation headrest offers precise adjustments for optimal head and neck support, enhancing patient comfort throughout any procedure. Additionally, it enhances dental access and operational efficiency, contributing to overall treatment success.



Closer to
Your Patients,
Closer to
Your Goals.



Belmont
CHAIRLOVE

**Are you prepared to elevate
your patient care experience
and strengthen connections
with your patients?**

Experience the Pro III dental chair — an embodiment of comfort, efficiency, and patient-centered care. Its innovative design brings you closer to your patients, supports your goals, and transforms every treatment room into a sanctuary of excellence.

Belmont CHAIRLOVE

Imagine a dental chair that not only meets the practical demands of dentistry but also enhances the patient experience at every touchpoint. The Pro III chair achieves this seamlessly, with features designed to foster deeper connections and improve patient comfort. From effortless patient access to seamless communication, every aspect of the Pro III is crafted to create a personalized and memorable interaction.

As your patients relax into the plush cushions and experience the smooth movement provided by the Pro III, they'll feel valued, understood, and well-cared for. And as a dentist, you'll appreciate the chair's user-friendly design, which enhances your ability to perform procedures efficiently without sacrificing patient comfort.

The Pro III isn't just a chair — it's a symbol of your commitment to excellence and your dedication to delivering the highest standard of care. It stands as a testament to your unwavering pursuit of success and fulfillment in your practice.

With the Pro III chair, you're investing in more than just equipment; you're fostering deeper connections with your patients and creating a brighter future for your practice. Get ready to elevate your practice with the Pro III dental chair — where every patient interaction leaves a lasting impact.



Belmont Quality.

Proven Reliability


Our hydraulic cylinders ensure smooth, precise movements essential for patient comfort and effective procedures. For decades, our designs have passed the rigorous 100,000-cycle durability test, assuring long-lasting performance.

By using fewer components and manufacturing in-house, we minimize malfunctions and ensure consistent quality, offering unparalleled precision and reinforcing Takara Belmont's status as a global leader in dental equipment.

Hydraulic Cylinder High Precision

Tolerance
Within 0.00083"

Tolerance
Within 0.00098"




**Customized for your workspace,
one stitch at a time.**



ULTRALEATHER®

Ultraleather offers luxury and durability, ideal for high-traffic areas. It maintains elegance and performance, even under demanding use.



ALLANTE FREE® STANDARD VINYL

Allante Free combines eco-friendly, PVC-free design with leatherlike finesse, offering a sustainable luxury option.



MAVERICK

Maverick vinyl blends durability and a vibrant color palette, enhancing any décor with a bright and modern aesthetic.



02

Evogue Delivery System

Evogue Cart Mobile Delivery System

**Experience unparalleled
versatility and convenience with
the Evogue cart.**

Meet the Evogue cart mobile-style delivery system, designed to elevate efficiency and flexibility in dental practices. With its compact, mobile design, it ensures effortless movement and positioning, optimizing space and workflow, while providing easy access to essential tools and instruments.





Streamlined Efficiency for Modern Practices.

Optimize Space & Workflow

The Evogue cart is an essential asset for small clinics and treatment rooms. Its compact design maximizes space, making it ideal for single-entry operatories where chair-mounted delivery isn't feasible. The cart's mobility keeps instruments within easy reach, enhancing workflow, supporting closer patient interactions, and fostering exceptional care in modern dental practices.



Enhanced Efficiency

The Evogue cart features height adjustability for easy access to handpieces and the instrument tray. Customizable to fit your needs, these carts can include built-in fiber optics, a micromotor, and a scaler, optimizing workspace efficiency. An optional stainless-steel tray provides a convenient location to place the rheostat, helping to keep your treatment room organized and clutter-free.





Evogue Cabinet-Mounted Delivery System

Evogue
Rear Delivery



The Evogue cabinet-mounted delivery system delivers space efficiency with two mounting options—rear and side delivery. Its discreet placement behind the patient helps reduce pre-treatment anxiety and improves the overall patient experience.

The Evogue cabinet-mounted delivery system acts as the operational hub of your treatment room, seamlessly integrating essential instruments and functions to enhance efficiency. Its discreet design reduces visual clutter, alleviating patient anxiety while optimizing your workflow for maximum accessibility.

Rear-Mounted Doctors & Assistants Delivery



Rear-Mounted Assistants Instruments



Side-Mounted Doctors Delivery



Compassion Meets Efficiency.

**How can your practice create
an environment of trust
and comfort for every patient?**

Dental care goes beyond procedures—it's about fostering trust, building relationships, and ensuring patients feel valued from the moment they step in.

The Evogue cabinet-mounted delivery system plays a pivotal role in shaping that experience, blending thoughtful design with seamless functionality to enhance patient comfort and practitioner efficiency.

Imagine a treatment space free from visual clutter, where sleek, integrated design creates a calming atmosphere that helps reduce patient anxiety. The Evogue system discreetly positions essential tools out of sight yet keeps them within effortless reach, ensuring your focus remains on your patients.

For practitioners, the streamlined rear delivery system transforms the treatment area into an efficient and organized workspace. With tools conveniently behind the chair, the design supports smooth task transitions, helping you provide care with confidence and ease.

By the end of each appointment, patients leave feeling reassured, valued, and confident in their care. The Evogue system is more than a tool—it's an investment in creating meaningful patient connections and elevating your practice. Let the Evogue system redefine what compassionate care looks like, one patient at a time.





Belmont Equipment

101 Belmont Drive, Somerset, NJ 08873
(800) 223-1192 | belmontequip.com

Scan for All Brochures

