



Belmont Dental Chair Collection

Q5500
Dental Chair

Q3300
Dental Chair

PRO III 047
Dental Chair

B50N
Dental Chair



Index

Belmont Dental Chair Collection

01 **Q5500** Dental Chair

02 **Q3300** Dental Chair

03 **Pro III** Dental Chair

04 **B50N** Dental Chair

01

—

Q5500
Dental Chair

**The Only Chair
Worthy of Center Stage.**



Exceptional Craftsmanship, Unsurpassed Quality

The Quolis Q5500 dental chair sets the stage for excellence through exceptional craftsmanship and quality. Its sleek, refined design enhances the experience for both patients and practitioners, blending intuitive features, user-friendly controls, and unmatched comfort. Take your practice to the center stage of excellence.



Elevate Care, Maximize Efficiency

Elevate care and maximize efficiency with the Q5500 chair, designed for superior patient outcomes and enhanced procedures. Combining smooth operation, illuminated controls, and a contoured backrest, it ensures optimal comfort and accessibility. By enriching the dental experience for both patients and practitioners, the Q5500 chair becomes an essential addition to any practice.



Intuitive Controls

The Q5500 chair is equipped with membrane touchpad controls below each armrest for easy adjustments. Each touchpad includes buttons for adjusting seat height and backrest positioning, programmable preset positions, safety lock-out to prevent unintentional seat motion, and a button to unlock the 120-degree seat rotation for precise control.



Oral Cavity Access

The dual-articulating headrest rotates forward or backward to improve access to the patient's oral cavity. This allows precise positioning for viewing and treating the maxillary and mandibular arches, reducing strain during both routine and complex procedures.





Superior Access

The Q5500 chair's slim, contoured backrest provides optimal access to the patient's oral cavity, allowing close positioning for precise, effortless procedures. Its design improves working posture and minimizes discomfort during both routine and complex treatments.



Smooth Operation, Seamless Movements

The Q5500 chair's shockless solenoids ensure gentle and seamless movements, enhancing patient comfort and practitioner control. This advanced technology provides a reliable, high-performance chair tailored to the dynamic needs of modern dental practices.



Reassuring Qualities

The Q5500 chair features user-friendly, rotating armrests that facilitate smooth patient entry and exit, reducing strain for both patients and staff. Designed for frequent use and stability, these armrests accommodate patients of all ages and mobility levels.



Thoughtfully designed with your needs in mind.

Crafted to meet the highest standards of modern dentistry, our chairs offer four distinct styles, each tailored to your unique clinical environment and comfort needs. Each style is thoughtfully designed to cater to the ergonomic and aesthetic needs of today's dental practices.



Plush Wide



Seamless Wide



Plush Narrow



Seamless Narrow





Belmont
CHAIRLOVE

**Better Patient Care
Starts with
a Better Patient Chair.**


**Do you aspire to redefine the
patient experience and lead your
dental practice in exceptional care?**

Belmont CHAIRLOVE

With the Q5500 dental chair, you can turn your vision into reality. Patients will feel at ease with the Q5500 chair's thin, narrow backrest and precise adjustments, while you'll appreciate its user-friendly design and efficiency-boosting features.

Throughout your day, the chair's dual-articulating headrest will prove invaluable, providing unmatched support during lengthy procedures and enhancing overall comfort. Patient feedback will highlight the chair's impact, with praise for the contoured backrest's spine support and the dual-articulating headrest's relief from neck strain.

As patients leave with smiles, it will become evident that the Q5500 is more than just a chair—it's an essential part of your practice. By day's end, the Q5500 will seamlessly integrate into your workflow, enhancing patient comfort and your efficiency. With its user-friendly design and unwavering focus on comfort and accessibility, the Q5500 chair plays a pivotal role in elevating your practice to a new level of excellence.



Belmont Quality.

Proven Reliability

Our hydraulic cylinders ensure smooth, precise movements essential for patient comfort and effective procedures. For decades, our designs have passed the rigorous 100,000-cycle durability test, assuring long-lasting performance.

By using fewer components and manufacturing in-house, we minimize malfunctions and ensure consistent quality, offering unparalleled precision and reinforcing Takara Belmont's status as a global leader in dental equipment.




**Customized for your workspace,
one stitch at a time.**



ULTRALEATHER®

Ultraleather offers luxury and durability, ideal for high-traffic areas. It maintains elegance and performance, even under demanding use.



ALLANTE FREE® STANDARD VINYL

Allante Free combines eco-friendly, PVC-free design with leatherlike finesse, offering a sustainable luxury option.



MAVERICK

Maverick vinyl blends durability and a vibrant color palette, enhancing any décor with a bright and modern aesthetic.



02

Q3300
Dental Chair

Better Patient Care Starts with a Better Patient Chair.

Optimal Comfort, Practitioner Efficiency

The Evogue Q3300 dental chair is more than just a chair—it's a partner in delivering superior care. Designed for optimal patient comfort and practitioner efficiency, it offers advanced ergonomics and intuitive controls in a sleek, modern design.



**Comfort made easy,
care made better.**

The Q3300 chair enhances patient comfort with a slim backrest for optimal access and relaxation. Intuitive controls, smooth motion, and effortless rotation streamline adjustments for practitioners, boosting efficiency and allowing greater focus on patient care.



Convenient Swivel Controls

Each chair features accessible seat rotation controls on the backrest, with push-buttons that automatically unlock the seat for effortless swiveling. This feature enables flexible patient orientations, enhancing the practitioner's view and access to the oral cavity.



Oral Cavity Access

The dual-articulating headrest mechanism allows the cushion to be rotated forward or backward for better viewing and access to the maxillary and mandibular arches. This precise adjustment enhances visibility, reduces physical fatigue, and minimizes interruptions for smoother, more accurate procedures.



Reassuring Qualities

The armrests rotate outward to facilitate smoother patient entry and exit, alleviating strain for both patients and dental staff. Built for frequent use and stability, these armrests accommodate patients of all ages and mobility levels.



Optimal Access

The Q3300 chair features a slim, contoured backrest that provides optimal access to the patient's oral cavity. Its streamlined design enables close positioning, minimizing discomfort and enhancing working posture during treatments.





**Thoughtfully designed
with your needs in mind.**

Crafted to meet the highest standards of modern dentistry, our chairs offer four distinct styles, each tailored to your unique clinical environment and comfort needs. Every style is thoughtfully designed to cater to the ergonomic and aesthetic requirements of today's dental practices.

Belmont
CHAIRLOVE

The Center of Their Success.

At the heart of every successful dental practice stands more than just a chair—it's a trusted partner in your journey to provide exceptional care.



Meet the Evogue Q3300 dental chair—designed to be more than just a piece of furniture. It's an integral part of your practice's success story. Serving as a steadfast companion, it embodies trust and dependability, standing as a symbol of your commitment to excellence.

Elevating the standards of dental chairs, the Q3300 masterfully integrates innovative design with exceptional functionality. Its synchronized reclining mechanism effortlessly adjusts to guide patients into optimal treatment positions, while its sleek, accessible design provides ease of reach from every angle. These smart features not only simplify the dentist's workflow but also ensure that every patient receives optimal care with minimal adjustments, redefining efficiency and comfort in modern dental practice.

For generations, Belmont has created dependable dental chairs that serve as the cornerstone of successful practices. Upholding unwavering standards of care is vital to cultivating a loyal patient base and a thriving practice. The Q3300 is more than just a chair—it embodies a commitment to a higher standard of care. With the Q3300 by your side, rest assured in your ability to provide your patients with the pinnacle of comfort and care, bolstering your practice for years to come.




**Customized for your workspace,
one stitch at a time.**



ULTRALEATHER®

Ultraleather offers luxury and durability, ideal for high-traffic areas. It maintains elegance and performance, even under demanding use.



ALLANTE FREE® STANDARD VINYL

Allante Free combines eco-friendly, PVC-free design with leatherlike finesse, offering a sustainable luxury option.



MAVERICK

Maverick vinyl blends durability and a vibrant color palette, enhancing any décor with a bright and modern aesthetic.



03

Pro III
Dental Chair

A Chair that Supports You from Every Angle.

**Experience the ultimate in dental
comfort and functionality with the
Pro III dental chair.**

Introducing the Pro III dental chair—a game-changer in innovation, functionality, and comfort. Tailored for both practitioners and patients, it merges advanced ergonomics with user-friendly features, setting a new benchmark in dental care. The Pro III offers superior support, ensuring comfort and confidence at every appointment.



Folding Legrest & Seat Rotation


The Pro III chair elevates the dental experience with features that set it apart from the rest. Its easy patient entry, space-saving efficiency, and unique design elements, such as swing-out armrests, automatic legrest extension, and a dual-articulating headrest, surpass standard offerings, providing unmatched comfort for all patients.





Upright Position

The extendable legrest and folding armrests facilitate easy patient entry before smoothly transitioning to a traditional dental setup for treatment.



Reclined Position

Built-in foot control levers reduce floor clutter around the chair and allow operators to transition seated patients to the desired treatment position.



Exceptional Capabilities, All-Around Versatility

The Pro III chair offers unparalleled flexibility with its 210-degree rotation range, allowing effortless positioning for easy patient entry, relaxed consultations, and seamless examinations from multiple angles.





10:00 Position



Doctors and staff can adjust the chair to create a welcoming environment and facilitate face-to-face consultations, reducing patient anxiety.



9:00 Position



Doctors can pivot patients to an optimal position, enhancing interactions and fostering a comfortable trusting environment.



6:00 Position



The unique folding legrest allows doctors to access patients seated upright directly from the front, promoting a natural, relaxed patient position for effective jaw alignment assessment and restorative procedures.

Reassuring Qualities

The Pro III chair is designed for optimal accessibility, with rotating armrests that ease patient entry and exit, reducing strain for both patients and staff. Built for durability and stability, the chair supports patients of all ages and mobility levels, with upright capabilities and rotating armrests providing easy, unobstructed access.



Oral Cavity Access

The chair's dual-articulation headrest offers precise adjustments for optimal head and neck support, enhancing patient comfort throughout any procedure. Additionally, it enhances dental access and operational efficiency, contributing to overall treatment success.



Closer to
Your Patients,
Closer to
Your Goals.



Belmont
CHAIRLOVE

**Are you prepared to elevate
your patient care experience
and strengthen connections
with your patients?**

Experience the Pro III dental chair — an embodiment of comfort, efficiency, and patient-centered care. Its innovative design brings you closer to your patients, supports your goals, and transforms every treatment room into a sanctuary of excellence.

Belmont CHAIRLOVE

Imagine a dental chair that not only meets the practical demands of dentistry but also enhances the patient experience at every touchpoint. The Pro III chair achieves this seamlessly, with features designed to foster deeper connections and improve patient comfort. From effortless patient access to seamless communication, every aspect of the Pro III is crafted to create a personalized and memorable interaction.

As your patients relax into the plush cushions and experience the smooth movement provided by the Pro III, they'll feel valued, understood, and well-cared for. And as a dentist, you'll appreciate the chair's user-friendly design, which enhances your ability to perform procedures efficiently without sacrificing patient comfort.

The Pro III isn't just a chair — it's a symbol of your commitment to excellence and your dedication to delivering the highest standard of care. It stands as a testament to your unwavering pursuit of success and fulfillment in your practice.

With the Pro III chair, you're investing in more than just equipment; you're fostering deeper connections with your patients and creating a brighter future for your practice. Get ready to elevate your practice with the Pro III dental chair — where every patient interaction leaves a lasting impact.




**Customized for your workspace,
one stitch at a time.**



ULTRALEATHER®

Ultraleather offers luxury and durability, ideal for high-traffic areas. It maintains elegance and performance, even under demanding use.



ALLANTE FREE® STANDARD VINYL

Allante Free combines eco-friendly, PVC-free design with leatherlike finesse, offering a sustainable luxury option.



MAVERICK

Maverick vinyl blends durability and a vibrant color palette, enhancing any décor with a bright and modern aesthetic.




04

B50N
Dental Chair

A Trusted Choice for Growing Clinics.

**Start on the right path to success
with quality chairs that support
your growth.**

Set your practice up for success with the X-Calibur V B50N dental chair. Its reliable electro-hydraulic system provides seamless motion and precise positioning, ensuring easy patient access without compromising comfort or safety. The perfect blend of durability and design, this chair is essential for any growing practice.



B50N
Dental Chair



Intuitive Controls

Features include a standard foot control and built-in membrane switches for precise patient seat positioning. The foot control allows hands-free operation, offering manual adjustments, two programmable presets, last position recall, and auto-return for optimal workflow efficiency. Built-in membrane switches, accessible from both the doctor and assistant sides, provide intuitive controls for seat height and backrest adjustments, enhancing ease of use and patient comfort.



Oral Cavity Access

The Belmont headrest's dual-articulation mechanism provides precise adjustments for optimal patient positioning, enhancing visibility and access. This reduces physical fatigue and improves overall efficiency while ensuring patient comfort during extended treatments.



Enhanced Accessibility

Reassuring qualities like outward-rotating armrests and a safety shut-off plate make the B50N chair a reliable choice for both doctors and patients. These accessibility features ensure easy entry and exit while preventing accidental activation, safeguarding the operator and the chair during every procedure.



Removable Sling Option

The removable sling upholstery is designed for easy cleaning of high-touch surfaces and effortless maintenance in busy clinics. This practical feature also allows for easy replacement as needed, ensuring a longlasting, hygienic solution.

For specialized care, such as pediatric dentistry, the sling can be removed to allow closer access to the patient.



Thoughtfully designed with your needs in mind.

Crafted to meet the highest standards of modern dentistry, our chairs offer four distinct styles, each tailored to your unique clinical environment and comfort needs. Every style is thoughtfully designed to cater to the ergonomic and aesthetic requirements of today's dental practices.



Plush
without Slings



Plush
with Slings



Standard
without Slings



Standard
with Slings


Belmont
CHAIRLOVE

Success Starts Here.

Recognizing that every decision impacts the growth of your dental practice, the B50N dental chair is expertly designed to be your ideal partner, checking all the right boxes for your evolving needs.

Start your day with confidence, knowing the B50N chair combines exceptional quality with great value. Its sleek design, 60° seat rotation, and electro-hydraulic system ensure effortless adjustments and optimal patient positioning, enhancing workflow efficiency. The dual-articulating headrest minimizes strain and maximizes comfort, while rotating arms and ETL-certified safety features add convenience and peace of mind.

With the B50N, you're not just investing in a chair; you're investing in a foundation for success. Its combination of reliability, durability, and user-friendly features makes it the perfect choice for any growing practice.




**Customized for your workspace,
one stitch at a time.**



ULTRALEATHER®

Ultraleather offers luxury and durability, ideal for high-traffic areas. It maintains elegance and performance, even under demanding use.



ALLANTE FREE® STANDARD VINYL

Allante Free combines eco-friendly, PVC-free design with leatherlike finesse, offering a sustainable luxury option.



MAVERICK

Maverick vinyl blends durability and a vibrant color palette, enhancing any décor with a bright and modern aesthetic.





Belmont Equipment

101 Belmont Drive, Somerset, NJ 08873
(800) 223-1192 | belmontequip.com

Scan for All Brochures

