



# BDU

Delivery System



# BDU Delivery System

A true example of reliability, flexibility, and durability.



**BDU570**  
Swing Mount



**BDU500**  
Over The Patient

## Easy Placement

Balance arm with pneumatic brake places the delivery head where it's needed and locks it in place.

## Operatory Footprint

Small delivery head takes up minimal space with a flat top to secure a stainless steel tray and non-slip pad to lay out instrumentation.

## Integrated Features

Pre-wired to support unit mounted touchpad and other add-on accessories.

## Versatility

Mounts to Belmont's B50N and 037N dental chairs, also works with Bel-Halo, Clesta LED and Clesta Halogen Lights.

## Water System

City/bottle water system with water bottle quick disconnect attached to PMU or Swing Arm.

Belmont Equipment  
101 Belmont Drive  
Somerset, NJ 08873  
(800) 223-1192  
[www.belmontequip.com](http://www.belmontequip.com)

Takara Company Canada Ltd.  
2076 South Sheridan Way  
Mississauga, ON  
L5J 2M4, Canada  
(800) 268-5351  
[dental.takarabelmont.ca](http://dental.takarabelmont.ca)

 **Belmont**  
Delivering Value in Quality

@BelmontEquipUS

