

Features & Performance	Bel-Halo	Bel-Nova	Clesta LED
Number of LEDs in Light Array	10	6	5
Light Pattern Shape	Rectangle	Oval	Rectangle
3-Axis Rotation Light Head	✓	✓	✓
Shadowless Illumination	✓	✓	✓
Correct Color	✓	✓	✓
Detachable Light Cover	-	✓	-
Rotating Patient Mirror	Option	Option	-
System & Controls			
Touchless Sensor for On/Off Activation	✓	✓	✓
Stepless Brightness Control	✓	✓	✓
Composite Safe Mode	✓	✓	✓
Blinking Composite Safe Mode Indicator Ring	✓	-	-
Manual Sensor On/Off Switch	✓	✓	✓
Styles & Materials			
LED Light Bulbs	✓	✓	✓
Mounting Options	Straight Post Mounted, Curved Post Mounted (For PMU), Swing Mounted, Ceiling Mounted, Track Mounted, Cabinet Mounted, Wall Mounted	Straight Post Mounted, Curved Post Mounted (For PMU), Swing Mounted, Ceiling Mounted, Track Mounted, Cabinet Mounted, Wall Mounted	Straight Post Mounted, Curved Post Mounted (For PMU), Swing Mounted, Ceiling Mounted, Track Mounted
Warranty			
5-Year Warranty	✓	✓	✓