


PRE-INSTALLATION INSTRUCTIONS MOUNTING TEMPLATE FOR WALL PLATE

Level Top and Side Prior Securing to Wall



mm
(inches)

MODEL 096,071A,DX-073,DX-068

1032 (40-5/8) FROM FLOOR

Belmont reserves the right to modify the specifications

EJPB03B0

<p>BELMONT EQUIPMENT CORP. 101 Belmont Drive Somerset, New Jersey 08873 U.S.A.</p>	<p>TAKARA CO. CANADA LTD. 2076 S.Sheridan Way Mississauga, Ont. L5J2M4 Canada</p>	<p>TAKARA BELMONT CORPORATION 1-1, 2-Chome Higashi-shinsaibashi, Chuo-ku Osaka Japan 542-0083</p>
---	--	--