



***ECO - SYS***  
Series Cabinetry

# ■ D-Series

## Operatory Cabinetry

With curved edges and premium options, the D-Series cabinetry offers concealed storage and enhances fluid workspace productivity.



### ECO1 Doctor's Console

Supplying the Doctor with storage and function.

- Large workspace with backsplash
- Stainless steel pulls on doors
- Stainless steel sink with chrome faucet and chrome soap dispenser
- Delivery system mounting capabilities
- Left- and right-side sink configurations available
- 24" under-sink storage area
- Waste receptacle with basket

### ECO2 Wall Mount Storage

Utilize vertical storage for a space-saving solutions.

- Stainless steel handles/pulls
- Dual glove/tissue dispenser on door
- Cup and towel dispenser
- Upper storage equipped with 4 adjustable shelves
- Can be pieced based on operatory needs (ECO2U or ECO2L)



### ECO3 Doctor's Console

Places storage and organized dispensers for everyday supplies at one's reach.

- Laminated work surface with drop-in stainless steel sink
- Utilizes 100-pound Unigrass drawer system
- Waste receptacle with basket
- Stainless steels handles/pulls
- Wash station includes: sink, manual faucet, soap dispenser, and hardware
- Optional barrier drawer system dispenses sanitation wipes, chair covers, light and syringe sleeves, tray covers, and cover-all roll dispenser
- Choose from left- or right-handed sink configuration





## ECO6 Rear Treatment Console

Adds an additional workspace within the operator with the flexibility of concealing for greater access in the office.

- Pull-out sliding countertop
- Delivery system mounting capability
- Task light and clock/timer
- Adjustable shelves on upper portions
- Easy accessible side drawers allow for additional storage space.
- Choose between left- or right-hand configuration
- Optional pass-thru on upper storage

## ECO8 Sink Unit

A one-stop for cleanliness and concealed disposal of waste.

- Standard laminated work surface with drop-in stainless steel sink
- Wash station includes: faucet, sink, and soap dispenser
- Stainless steel handles/pulls on doors
- Waste receptacle with basket



## ECO9 Mobile Cart

Positions a workspace and storage at arm's reach.

- Choose worksurface slide-out orientation (front-back, left-right)
- Additional 9" storage bin beneath worksurface with electrical access grommet
- Stores neatly in storage bay of the ECO1 Doctor's Console
- Stainless steel handles/pulls on drawers
- 4 swivel casters for mobility and stability

All ECO-SYS cabinets come standard with Laminate Color option and Wilsonart® HD worksurface colors. Additional worksurface options are available as an upgrade.

### Standard Laminate



### Wilsonart® H.D. Worksurface



### LG Hausys HI-MACS® Worksurface (Upgrade)



Belmont Equipment  
101 Belmont Drive  
Somerset, NJ 08873  
(800) 223-1192  
www.belmontequip.com

Koken Manufacturing Company, Inc.  
A Division of Takara Belmont USA, Inc.  
1631 Martin Luther King Drive  
St. Louis, MO 63106  
T: (800) 325-7373  
F: (314) 446-7088

Takara Company Canada Ltd.  
2076 South Sheridan Way  
Mississauga, ON  
L5J 2M4, Canada  
(800) 268-5351  
dental.takarabelmont.ca