



IMPROVE YOUR PROCEDURAL SET UP
WITH THE EVOGUE FLEX

Ergonomics, durability and aesthetics combined together to create one of the most efficient and easy to use delivery systems available.



VERSATILE

Evogue Flex has a pivoting doctor control which can be positioned for side delivery, as shown, or directly over the patient's chest, for traditional continental-style delivery.

ERGONOMIC

Handpiece actuation arms utilize a unique counter balance adjustment mechanism which reduce strain and enhances operator comfort during patient treatment.

DURABLE


Quality materials and rugged construction help to ensure that each unit functions reliably and looks great for many years to come.

LONGER REACH. BETTER RESULTS.


Evogue Flex delivery systems route handpiece hoses externally from actuator arms, and through rotating tip pulley systems. This efficient design provides longer reach

of each handpiece, so the need for operator repositioning and/or delivery head repositioning are minimized during patient treatment.

EVOGUE Flex



Competitor



CENTRALIZE YOUR PROCEDURES

The Evogue Flex delivery head combines unique handpiece arm actuation systems and a pivoting doctor table, both of which result in greater accessibility and comfort for operators during patient treatment.



IMPROVED ACCESS, LESS STRAIN

Handpieces are activated by pulling hoses directly forward or from a side angle. This enables operators to minimize upper body motion and maintain their focus where it belongs – on the patient.



THE COLLECTION THAT YOU NEED



The Evogue Delivery System complements Belmont's robust product line of dental chairs, delivery systems, lights, x-ray and cabinets. Outfit your operatory with Belmont's signature quality products today.

Clesta LED Dental Light



Bel-Halo Dental Light



Flex Rear Cabinet



Quolis® Q5000 Dental Chair



Bel50N Dental Chair



Flex Side Cabinet



PHOT-xlIs LCD X-Ray



80B Dental Stool



81B Dental Stool



Belmont Equipment
101 Belmont Drive
Somerset, NJ 08873
(800) 223-1192
www.belmontequip.com

 **Belmont**
Delivering Value in Quality

@BelmontEquipUS

