



## Sterilization Centers

# ■ E 155

## Sterilization Center

Luxurious style. Quality performance.




### Hands-free Operations

Color-coded pneumatic switches to operate hands-free waste drawer (green), unsterile upper storage doors (red), and hands-free tower doors (blue).



### Sterile

Red light storage for soiled instrument cassettes and sanitized instrument storage with blue light to indicate items have been processed.



### Safe Environment

Pneumatic hands-free waste drawer with sharps easily manages waste materials based on safety guidelines. The push-out opening mechanism allows the user to open and close the drawer without any cross-contamination.



### Accessibility

Tapmaster™ foot control can be found underneath the sink for hands-free operations.

# ■ D96

## Sterilization Centers

Incorporates a modular design with customizing components to deliver simplistic function.



### Sink Usability

Tapmaster™ kick pedal can be found underneath the sink for hands-free operations



### Workspace Components

Standard laminate countertop with drop-in stainless steel sink and chrome faucet with soap dispenser.



### Steri-Tower

Sterilization tower utilizes wrap drawers with slide-out shelves to keep tools accessible and within reach.



### Clean Atmosphere

Pull-out trash module keeps waste out of sight for a cleaner look.



Belmont Equipment  
101 Belmont Drive  
Somerset, NJ 08873  
(800) 223-1192  
[www.belmontequip.com](http://www.belmontequip.com)

Koken Manufacturing Company, Inc.  
A Division of Takara Belmont USA, Inc.  
1631 Martin Luther King Drive  
St. Louis, MO 63106  
T: (800) 325-7373  
F: (314) 446-7088

Takara Company Canada Ltd.  
2076 South Sheridan Way  
Mississauga, ON  
L5J 2M4, Canada  
(800) 268-5351  
[dental.takarabelmont.ca](http://dental.takarabelmont.ca)